

From X-rays to Point Groups, one step at a time

CCCW25

June 20, 2025

Andreas Decken (via Brian Patrick via Kate
Markzenko via Mike Katz)

Designer

Sorry, no design ideas for this slide.


When we have design ideas, we'll show them to you right here.

From X-rays to Point Groups, one step at a time

CCCW25

June 20, 2025


Andreas Decken (via Brian Patrick via Kate
Markzenko via Mike Katz)




Introduction

Understanding the Structure of Materials

- We use **crystallography** to understand the structure of **crystalline** solids.
 - **Crystallography:** The study of *crystal structures*.
 - **Crystal structure:** A description of the ordered arrangement of atoms, ions, or molecules in a crystalline solid.
 - **Crystalline solid:** A material whose constituents are arranged in a high order, forming a crystal lattice
 - **Crystal lattice:** 3D arrangements of atoms in a crystal
 - **Crystal:** or crystalline solid...
 - Why do you care about the structure of your crystalline solids?



Electromagnetic Radiation




- Observing small objects requires short wavelength
- X-rays can penetrate matter more easily than visible light.

A New Kind of Rays


■ Wilhelm Conrad Röntgen

He gave his first public lecture on X-rays in January 1896 and showed the rays' ability to photograph the bones within living flesh. A few weeks later in Canada, an X-ray was used to find a bullet in a patient's leg.

Röntgen deliberately didn't patent his discovery, feeling that scientific advances belonged to the world and should not be for profit.




X-Ray Tube



- Generation of lots of heat
- More powerful variants are Rotating X-ray sources
- Liquid Metal Stream
- Microfocus Tubes

Source: Andreas Eckmüller


Generation of Characteristic Radiation



- Incoming electron knocks out an electron from the inner shell of an atom.
- Designation K,L,M correspond to shells with a different principal quantum number.

Emission Spectrum of an X-Ray Tube: Close-up of $K\alpha$

Bremsstrahlung
radiation



Characteristic Wavelength

The X-rays we use for diffraction have a wavelength of 0.7-1.5 Å.


- Cu radiation (1.54 Å)
 - Great for organic compounds and chiral molecules.
 - Diffracted spots are more widely spread; great for large unit cells.
- Mo radiation (0.71 Å)
 - Great for inorganic compounds.
 - Fast data collection.

Element	Wavelength (λ) Å ($K\alpha_1$)
Ag	0.5594
Mo	0.7097
Cu	1.5405
Co	1.7889
Fe	1.9360
Cr	2.2896

PXRD Analysis

X-Ray Scattering

How does electromagnetic radiation interact with matter?



- X-rays interact with the electrons in an atom, causing them to change directions and scatter.
 - **Elastic:** Same wavelength and energy.
 - Inelastic: Different wavelength and energy.
- With multielectron atoms, the whole electron cloud oscillates.

X-Ray Scattering


What happens when we have more than one atom?



The scattered waves from each atom can **interfere** with each other.

In Phase


Wave add together



Constructive Interference

180° Out of Phase

Wave cancel each other




Destructive Interference

Source:



X-Ray Scattering

Young's Double Slit Experiment




- Light passing through the slits creates two separate wavefronts that can interfere either:
 - Constructively (amplitudes of waves add).
 - Destructively (amplitudes of waves cancel).


Bragg's Law

Model

- Crystal is a set of discrete parallel planes.





Derivation of Bragg's Law




$$\sin \Theta = n \lambda / 2 d$$

$$n \lambda = 2d \sin \Theta$$

Source: Patrick Woodward





Source: Taylor & Lipson, Optical Transforms








••


Source: Taylor & Lipson, Optical Transforms




Source: Taylor & Lipson, Optical Transforms



Source: Taylor & Lipson, Optical Transforms



Planes and Distances




- We can draw an infinite number of parallel planes through a crystal.
- The family of planes are a distance d apart from each other.
- Only when we satisfy Bragg's Law, do we see "diffraction" from the planes.

$$2d \cdot \sin(\theta) = n \cdot \lambda$$

Planes and Distances

Miller Planes


- In 1839, William Hallowes Miller came up with a notation for these planes.
- We don't need to consider the whole crystal. We just need to look at the unit cell.



Planes and Distances

Miller Planes

- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a , b , c).




Planes and Distances

Miller Planes

- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a , b , c).

The orange plane:

- Intersects **a** at 0.5
- Intersects **b** at 0.5
- Intersects **c** at 1



Planes and Distances

Miller Planes


- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a , b , c).

The orange plane:

- Intersects \mathbf{a} at 0.5
- Intersects \mathbf{b} at 0.5
- Intersects \mathbf{c} at 1

We refer to these planes as $(h k l)$.

- h is the inverse of the fractional coordinate along \mathbf{a} .
- k is the inverse of the fractional coordinate along \mathbf{b} .
- l is the inverse of the fractional coordinate along \mathbf{c} .




Planes and Distances

Miller Planes

- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a , b , c).

The orange plane:

- Intersects \mathbf{a} at 0.5
 - Intersects \mathbf{b} at 0.5
 - Intersects \mathbf{c} at 1
- h is the inverse of the fractional coordinate along \mathbf{a} .
 - $h = 1/0.5 = 2$
 - k is the inverse of the fractional coordinate along \mathbf{b} .
 - $k = 1/0.5 = 2$
 - l is the inverse of the fractional coordinate along \mathbf{c} .
 - $l = 1/1 = 1$




Planes and Distances

Miller Planes

- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a , b , c).

The orange plane:

- Intersects \mathbf{a} at 0.5
 - Intersects \mathbf{b} at 0.5
 - Intersects \mathbf{c} at 1
- h is the inverse of the fractional coordinate along \mathbf{a} .
 - $h = 1/0.5 = 2$
 - k is the inverse of the fractional coordinate along \mathbf{b} .
 - $k = 1/0.5 = 2$
 - l is the inverse of the fractional coordinate along \mathbf{c} .
 - $l = 1/1 = 1$




Planes and Distances

Miller Planes

- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a, b, c).

The purple plane:

- Intersects a at 1
- Intersects b at 0.5
- Intersects c at 0.33




Planes and Distances


Miller Planes


- Every plane that we can draw in the unit cell, intersects at some fraction of the way along the unit cell axis (a, b, c).


The purple plane:


- Intersects a at 1
- Intersects b at 0.5
- Intersects c at 0.33






















Practice




Practice




(200)




Practice




(200)




(111)




Practice



(200)



(111)



(044)

Miller Plane Calculator, Michael Katz

MEMORIAL UNIVERSITY

Home About Us Research Publications Teaching Media Contact

Dr. Michael Katz - Porous Materials

Miller Planes

Below you will find a Miller plane viewer. Miller planes can be viewed as planes that slice through a unit cell. It can be used to visualize crystal growth directions. It can also be used to consider how X-ray "diffraction" (actually scattering) occurs off of a plane. We define each plane using a set of indices, h , k , and l . Each index represents the inverse of how far along a unit cell the plane passes through along x , y , and z respectively. For example, if $h = 1$, then the plane intersects the x -axis at a . Furthermore, if $h = 2$, then the plane intersects the x -axis at half of a ($t/h = 1/2$). We label planes as " $(h k l)$ ". So the plane $(1 1 0)$ intersects the x -axis at a , the y -axis at b , and it runs parallel to the z -axis $t/l = 1/0 = \infty$.

Use the viewer below to explore the viewer. Be careful with negative values of h , k , and l .

Miller Plane 1 (green):

$h:$

$k:$

$l:$

Miller Plane 2 (purple):

Enable

$h:$

$k:$

$l:$

a length:

b length:


c length:

Apply Changes

Contact Details

Department of Chemistry,
Memorial University,
St. John's, NL
Canada


Phone: 1-709-864-8745
Email: mkatz@mun.ca
Website: www.KatzResearchGroup.com



If we know the unit cell, and how it is oriented on the instrument, we can calculate exactly where a reflection should strike the detector!

Crystal Systems and Crystallographic Point Groups

Andreas Decken via Paul D. Boyle



Crystal Systems and Crystallographic Point Groups

Andreas Decken via Paul D. Boyle



Crystal Systems and Crystallographic Point Groups

Andreas Decken via Paul D. Boyle











1820 × 1212





Your Data is Precious

How many data points do we need to solve the following equation?

$$y = ax + b$$

$$y = ax^2 + bx + c$$


In an X-ray diffraction experiment we obtain a limited amount of data

If we want precise bond length/angles we want to optimize
how we use our data




Source:
NDK Synth.
Quartz


Crystal Systems




Edges and angles




Cubic
 $a = b = c$
 $\alpha = \beta = \gamma = 90^\circ$




Tetragonal
 $a = b \neq c$
 $\alpha = \beta = \gamma = 90^\circ$




Orthorhombic
 $a \neq b \neq c$
 $\alpha = \beta = \gamma = 90^\circ$




Monoclinic
 $a \neq b \neq c$
 $\alpha = \gamma = 90^\circ \neq \beta$



Hexagonal
 $a = b \neq c$
 $\alpha = \beta = 90^\circ, \gamma = 120^\circ$



Rhombohedral
 $a = b = c$
 $\alpha = \beta = \gamma \neq 90^\circ$



Triclinic
 $a \neq b \neq c$
 $\alpha \neq \beta \neq \gamma \neq 90^\circ$

Crystal Systems

Crystal System	Metric Constraints
Triclinic	$a \neq b \neq c; \alpha \neq \beta \neq \gamma$
Monoclinic	$a \neq b \neq c; \alpha = \gamma = 90^\circ, \beta \neq 90^\circ$
Orthorhombic	$a \neq b \neq c; \alpha = \beta = \gamma = 90^\circ$
Tetragonal	$a = b \neq c; \alpha = \beta = \gamma = 90^\circ$
Trigonal rhombohedral setting hexagonal setting	$a = b = c; \alpha = \beta = \gamma \neq 90^\circ$ $a = b \neq c; \alpha = \beta = 90^\circ, \gamma = 120^\circ$ $a = b \neq c; \alpha = \beta = 90^\circ, \gamma = 120^\circ$
Hexagonal	$a = b \neq c; \alpha = \beta = 90^\circ, \gamma = 120^\circ$
Cubic	$a = b = c; \alpha = \beta = \gamma = 90^\circ$

Metric Constraints

NOTE: “=” means “constrained to be equal”

“≠” means “not constrained to be equal to”
it does not mean “not equal to”

Systems – Point Groups – Space Groups

System

Triclinic

Monoclinic

Orthorhombic

Tetragonal

Trigonal/rhombohedral

Hexagonal

Cubic



Crystallographic Point Groups

A group of symmetry operations all for which leave at least one point unmoved. It maps a point lattice onto itself.

In 3D, the symmetry elements are restricted to:


- Proper rotation axes (n)
- Mirror planes (m)
- Inversion centre ($\bar{1}$, or no explicit symbol)
- Rotary inversion axes (\bar{n})

Proper Rotation Axes


- Rotation about an axis by $360^\circ/n$
- Doesn't change handedness of object




2 fold
'Diad'



3 fold
'Triad'



4 fold
'Tetrad'




6 fold
'Hexad'

- Only n-fold axes where $n = 1, 2, 3, 4, 6$ are allowed for space filling 3 dimensional objects


Mirror plane

- Creates a reflected object
- Changes handedness of object





Inversion Centre

- Transforms x, y, z into $\bar{x}, \bar{y}, \bar{z}$
- Changes handedness of object



No cows, no hands


Proper Rotation Mirror



Source: Stout & Jensen

Rotary Inversion Axes

- Rotation of $360^\circ/n$ followed by inversion
- Changes handedness of object
- $\bar{1}$ is equivalent to an inversion centre
- $\bar{2}$ is equivalent to a mirror plane



Source: Stout & Jensen

Crystallographic Point Groups

32 unique crystallographic point groups are obtained from combining the various allowed rotation axes, mirror planes, and inversions

11 of the 32 crystallographic point groups are ***centrosymmetric*** and are the ***Laue Groups***

Laue Groups describe the symmetry of the diffraction pattern as determined from the observed intensities

Determining the Laue Group is an important step when solving a crystal structure

Systems – Point Groups – Space Groups

System	Point Group
Triclinic	$\bar{1}$ $\bar{1}$
Monoclinic	2 m $2/m$
Orthorhombic	222 $mm2$ mmm
Tetragonal	4 $\bar{4}$ $4/m$ 422 $4mm$ $\bar{4}2m$

4/mmm

Trigonal/rhombohedral	$\frac{3}{\bar{3}}$
	32
	$\frac{3m}{\bar{3}m}$
Hexagonal	$\frac{6}{\bar{6}}$
	<i>6/m</i>
	622
	<i>6mm</i>
	$\bar{6}m\bar{2}$
	<i>6/mmm</i>
Cubic	23
	<i>m3</i>
	432
	$\bar{4}3m$
	<i>m3m</i>


Hermann-Mauguin Symmetry Notation

- Spectroscopists use Schoenflies notation to describe symmetry (e.g. C_{2v} , D_{4h})
- Crystallographers use Hermann-Mauguin notation (International notation)
- Was introduced by Carl Hermann in 1928
- Modified by Charles-Victor Mauguin in 1931
- Adopted for the 1935 edition of the *International Tables for Crystallography*
- Hermann-Mauguin notation is preferred for crystallography
 - Easier to add translational symmetry elements
 - Directions of symmetry axes are specified


Crystal family	Crystal system	Hermann-Mauguin		Shubnikov ^[2]	Schoenflies	Orbifold	Coxeter	Order
		(full)	(short)					
Triclinic		1	1	1	C_1	11	$[]^+$	1
		$\bar{1}$	$\bar{1}$	$\tilde{2}$	$C_i = S_2$	\times	$[2^+, 2^+]$	2
Monoclinic		2	2	2	C_2	22	$[2]^+$	2
		m	m	m	$C_s = C_{1h}$	*	$[]$	2
		$\frac{2}{m}$	2/m	$2 : m$	C_{2h}	2^*	$[2, 2^+]$	4
Orthorhombic		222	222	$2 : 2$	$D_2 = V$	222	$[2, 2]^+$	4
		mm2	mm2	$2 \cdot m$	C_{2v}	*22	$[2]$	4
		$\frac{2}{m} \frac{2}{m} \frac{2}{m}$	mmm	$m \cdot 2 : m$	$D_{2h} = V_h$	*222	$[2, 2]$	8
Tetragonal		4	4	4	C_4	44	$[4]^+$	4
		$\bar{4}$	$\bar{4}$	$\tilde{4}$	S_4	$2\times$	$[2^+, 4^+]$	4
		$\frac{4}{m}$	4/m	$4 : m$	C_{4h}	4^*	$[2, 4^+]$	8
		422	422	$4 : 2$	D_4	422	$[4, 2]^+$	8
		4mm	4mm	$4 \cdot m$	C_{4v}	*44	$[4]$	8
		$\bar{4}2m$	$\bar{4}2m$	$\tilde{4} \cdot m$	$D_{2d} = V_d$	2^*2	$[2^+, 4]$	8
		$\frac{4}{m} \frac{2}{m} \frac{2}{m}$	4/mmm	$m \cdot 4 : m$	D_{4h}	*422	$[4, 2]$	16
		3	3	3	C_3	33	$[3]^+$	3
Trigonal		$\bar{3}$	$\bar{3}$	$\tilde{6}$	$C_{3i} = S_6$	$3\times$	$[2^+, 6^+]$	6
		32	32	$3 : 2$	D_3	322	$[3, 2]^+$	6
		3m	3m	$3 \cdot m$	C_{3v}	*33	$[3]$	6
		$\bar{3} \frac{2}{m}$	$\bar{3}m$	$\tilde{6} \cdot m$	D_{3d}	2^*3	$[2^+, 6]$	12
		6	6	6	C_6	66	$[6]^+$	6
Hexagonal		$\bar{6}$	$\bar{6}$	$3 : m$	C_{3h}	3^*	$[2, 3^+]$	6
		$\frac{6}{m}$	6/m	$6 : m$	C_{6h}	6^*	$[2, 6^+]$	12
		622	622	$6 : 2$	D_6	622	$[6, 2]^+$	12
		6mm	6mm	$6 \cdot m$	C_{6v}	*66	$[6]$	12
		$\bar{6}m2$	$\bar{6}m2$	$m \cdot 3 : m$	D_{3h}	*322	$[3, 2]$	12
		$\frac{6}{m} \frac{2}{m} \frac{2}{m}$	6/mmm	$m \cdot 6 : m$	D_{6h}	*622	$[6, 2]$	24
		23	23	$3/2$	T	332	$[3, 3]^+$	12
		$\frac{2}{m} \bar{3}$	$m\bar{3}$	$\tilde{6}/2$	T_h	3^*2	$[3^+, 4]$	24
Cubic		432	432	$3/4$	O	432	$[4, 3]^+$	24
		$\bar{4}3m$	$\bar{4}3m$	$3/\bar{4}$	T_d	*332	$[3, 3]$	24
		$\frac{4}{m} \bar{3} \frac{2}{m}$	$m\bar{3}m$	$\tilde{6}/4$	O_h	*432	$[4, 3]$	48



TRICLINIC
Laue group $\bar{1}$




$1 (C_1)$




$\bar{1} (C_i)$


MONOCLINIC
Laue group $2/m$



$2 (C_2)$




$m (C_s)$




$2/m (C_{2h})$


ORTHORHOMBIC
Laue group mmm



$222 (D_2)$




$mm2 (C_{2v})$




$4 (C_4)$


TETRAGONAL
Laue group $4/m$



$4 (C_4)$




$\bar{4} (S_4)$




$4/m (C_{4h})$


TETRAGONAL
Laue group $\frac{4mm}{m}$




$422 (D_4)$




$4mm$ (C_{4v})




$\bar{4}2m$ (D_{2d})




$\frac{4}{m}mm$ (D_{4h})




6 (C_6)




622 (D_6)




$\bar{6}$ (C_{3h})




$6mm$ (C_{6v})



$6/m$ (C_{6h})




$6m2$ (D_{3h})




$\frac{6}{m}mm$ (D_{6h})

TRIGONAL Laue group $\bar{3}$




3 (C_3)




$\bar{3}$ (S_6)


TRIGONAL
Laue group $\bar{3}m$



$32 (D_3)$




$3m (C_{3v})$




$\bar{3}m (D_{3d})$

CUBIC
Laue group $m3$




$23 (T)$




$m3 (T_h)$


CUBIC Laue group $m3m$



$432 (O)$



$\bar{4}3m (T_d)$



$m3m (O_h)$

Source: Ladd & Palmer

Choosing the Correct Crystal System

- Do not assume the metric relations indicate the correct point group and crystal system!!!
- Correctly identify the Laue group symmetry of the diffraction pattern (equivalent intensities, R_{sym})
- The Laue symmetry indicates the crystal system of your sample
- Correct Laue group assignment narrows space group choices

Space Groups

- Space groups vs Point groups
 - Point groups describe symmetry of isolated objects
 - Space groups describe symmetry of infinitely repeating space filling objects
- Space groups include point symmetry elements

Space Groups

- Space groups vs Point groups
 - Point groups describe symmetry of isolated objects
 - Space groups describe symmetry of infinitely repeating space filling objects
- Space groups include point symmetry elements
- After the break - yeah

